

# Yeast Extract Sample Product Data Sheet (Catalog No. 12110-SM)

## Information for Ordering Product

Catalog No.	Quantity
12110-05	500 g
12110-1	1 kg
12110-5	5 kg
12110-10	10 kg

**This product is for research purposes only. Not for diagnostic use.**

### Description

Sample is provided as 10 g of powdered Yeast Extract.

Yeast Extract is the water soluble portion of autolyzed yeast with intact B-complex vitamins. It is ideal for general bacteriological use with a variety of microorganisms. Yeast Extract is guaranteed BSE compliant.

---

### Storage Conditions

Store at 20-25 °C for up to 1 year.

---

### Other Growth Media Products available from Mo Bio Laboratories, Inc.

<u>Product Description</u>	<u>Catalog No.</u>
TB DRY™ Powder Growth Media (1 kg)	12105-1
LB Broth Powder Growth Media (1 kg)	12106-1
LB Agar Powder Growth Media (1 kg)	12107-1
LB Broth (Lennox) Powder Growth Media (1 kg)	12108-1
LB Agar (Lennox) Powder Growth Media (1 kg)	12109-1
Soybean-Casein Digest Medium (TSB), USP (1 kg)	12114-1
Soybean-Casein Digest Agar Medium (TSA), USP (1 kg)	12115-1
Tryptone (1 kg)	12111-1
Agar (1 kg)	12112-1

---

### Ordering Information

If you received this product from one of our distributors, please call them again to order.  
For the distributor nearest you, please go to our web site at [www.mobio.com/distributors](http://www.mobio.com/distributors)

Manufactured by Mo Bio Laboratories, Inc. 2746 Loker Avenue West, Carlsbad, CA 92008  
Ph: 760-929-9911 800-606-6246 Fax: 760-929-0109 E-mail: [orders@mobio.com](mailto:orders@mobio.com)

Version 03212005

**Technical Support: 1-800-606-6246 or 1-760-929-9911 E-mail: [technical@mobio.com](mailto:technical@mobio.com)**